

**Backed  
by the  
Best**

## Stay Moving with Regular Maintenance



### A Mobility Expert's Guide to Maximizing the Life and Trade-In Value of Your Wheelchair Accessible Vehicle

A new mobility vehicle is more than just a big investment. It represents the independence to get where you need to go and the freedom to choose where you want to be each day.

That's why protecting your BraunAbility product with regular maintenance, both professionally, and at home, is so critical. With a little care, your accessible vehicle will be your gateway to the road for years to come while preserving maximum trade-in value. You'll have peace of mind knowing your investment and independence are well protected.

#### Professional Maintenance Schedule

We've created this easy checklist so you can mark your calendar in advance.

##### 5,000 Miles Or 6-Month Intervals

Perform specific inspection, maintenance and lubrication procedures.

###### 6-Month Dealer Maintenance Checklist

- ☐ Inspect ramp access door system operation and alignment; clean and lube system
- ☐ Inspect ramp system operation; clean and lube system
- ☐ Inspect kneel system operation; clean and lube system (if applicable)
- ☐ Inspect electronic control system for current software (if applicable)
- ☐ Inspect undercarriage and ground effects
- ☐ Inspect and maintain wheelchair tiedowns and occupant restraints

##### 30,000 Miles Or 3-Year Intervals

Perform specific inspection, maintenance and lubrication procedures.

###### 3-Year Dealer Maintenance Checklist

- ☐ Door rollers (if applicable); Inspect and replace as necessary
- ☐ Kneel chain (if applicable); Inspect and replace as necessary
- ☐ Foldout ramp chain (if applicable); Inspect and replace as necessary
- ☐ \*Infloor ramp belt and rollers (if applicable); Inspect and replace as necessary

\* Infloor ramp motor and drive components are accessible only by removing the top cover. Refer to applicable service manual.

*Continued*

# Stay Moving with Regular Maintenance



## Choosing The Right Professional

You purchased a BraunAbility mobility vehicle because it was specially designed for your needs. The same care should go into choosing the right professional to help you maintain it.

Your local BraunAbility service technician is specially trained to understand not just the vehicle, but the conversion components and how they are designed to operate. To find your closest Mobility Service Excellence (MSE) trained technician, call our Customer Care line at 1-800-488-0359.

## Protecting Your Investment at Home

Regular cleaning and inspections at home between maintenance visits will ensure your van is running strong when you need it most. Follow the guide below:

## At-Home Maintenance Quick Reference Chart

### Vehicle Parts

**Start Vehicle:** listen for stutters or hesitations

**Run Vehicle:** 20-30 minutes twice a week

**Battery:** Open hood, examine battery for debris on cables & contact points

**Wipers:** Check/Replace

**Tires:** check pressure (found in manual or inside car door tag) and inspect tread wear

### Mobility Parts

**Lower Door Track:** check for/remove debris\*

**Ramp:** Extend, lubricate\*\* hinge points, tighten bolts

**Door Rollers:** Check for/remove debris\*, lubricate\*\*

**Tie-Down Tracks:** Test, if slow to adjust/move: check for/remove debris\*

\*Debris and dirt can be dislodged with a brush or swept up with a vacuum hose.

\*\*Use light lubricant or dry lubricant spray.

## Know the Back-ups

In the rare instance that you experience a part failure, backup systems are in place to ensure your independence is protected. Your local mobility dealer can walk you through manual operations in advance, so you know exactly what to do in a pinch.